

Determining the Phosphate Content in Fertilizers

Introduction:

Fertilizers are widely used to improve the growth of crops. Many of these fertilizers contain phosphates.. This experiment will determine the amount of phosphorous present as phosphate.

Purpose:

The purpose of this experiment is to determine the amount of phosphorous present in a soluble fertilizer in the form of the phosphate ion.

Equipment/ Materials

HACH kit and supplies for the Phosphorus, Reactive (Orthophosphate) test (number 490)

10 mL sample cells

PhosVer3 phosphate Powder Pillow

Standard phosphate solution $\text{Na}_3\text{PO}_4 \cdot 12\text{H}_2\text{O}$ 50.0 mg/L

Volumetric flasks (25 mL)

Pipet (10 mL)

Pipet bulb

Kim wipes

Automatic pipet adjustable

Solutions of fertilizers containing phosphate

Safety:

Wear an apron and goggles

Procedure:

1. Use the 50 mg/L Standard solution to prepare 25 mL of standards containing 2.4 mg/L, 1.0 mg/L, 0.75 mg/L and 0.5 mg/L of phosphate. Also prepare a blank. The blank in this case is a 10 mL sample cell of distilled water.
2. Measure the concentration of each of these standards using the DR 2400 from the HACH kit.
 - a. Touch favorite programs
 - b. Select program 490 P React. PV
 - c. Touch Start
 - d. Fill a round 10 mL sample cell with 10 mL of a standard
 - e. Add the contents of one PhosVer 3 phosphate Powder Pillow to the standard. Immediately cap and invert to mix. This becomes the prepared sample
 - f. Touch the clock in the lower left side of the screen. This is the timer icon.
 - g. Touch Start. The clock begins to count down as a two-minute reaction period is started. Place the prepared sample vial on the lab table.

- h. Wipe the blank and place it into the cell holder.
 - i. When the timer beeps, Touch Zero. The display will show 0.0 mg/L PO₄³⁻.
 - j. Within one minute after the timer beeps, (indicating completion of the Zero reading), wipe the prepared sample and place it into the cell holder.
 - k. Touch Read.
 - l. Results will appear in mg/L PO₄³⁻. Record data
 - m. Repeat until all standards have been measured. Record data.
3. Obtain a sample of fertilizer solution
 4. Repeat step 2 using the fertilizer solution. Record data
 5. Repeat step 2 using additional fertilizers. Record data

Data:

Standards

Prepared Concentration (mg/L)	HACH meter reading (mg/L)
2.4	
1.0	
0.75	
0.5	

Fertilizer

Kind of fertilizer used _____

HACH meter reading _____

Kind of fertilizer used _____

HACH meter reading _____

Kind of fertilizer used _____

HACH meter reading _____

Calculations:

1. Using the equation $C_1V_1 = C_2V_2$, calculate the volume of 50.0 mg/L PO₄³⁻ solution needed to prepare 25 mL of a solution containing 2.4 mg/L PO₄³⁻

2. Using the equation $C_1V_1 = C_2V_2$, calculate the volume of 50.0 mg/L PO_4^{3-} solution needed to prepare 25 mL of a solution containing 1.0mg/L PO_4^{3-}

3. Using the equation $C_1V_1 = C_2V_2$, calculate the volume of 50.0 mg/L N solution needed to prepare 25 mL of a solution containing 0.75 mg/L PO_4^{3-}

4. Using the equation $C_1V_1 = C_2V_2$, calculate the volume of 50.0 mg/L PO_4^{3-} solution needed to prepare 25 mL of a solution containing 0.50 mg/L PO_4^{3-}

5. Graph your data. Plot the prepared concentrations on the x axis and the HACH meter readings on the Y axis

6. Using the calibration curve constructed above and the HACH meter reading, calculate the concentration of PO_4^{3-} in your fertilizer.

Questions

1. Why was a calibration curve constructed.?
2. List three errors that would cause phosphate readings to be too low.
3. List two errors that would cause phosphate readings to be too high.